



THROMBO TIMER Plus

Activated Clot Detection System

- ⌚ **THROMBO TIMER PLUS** is an Activated Clot Detection Machine designed to perform tests like PT, APTT.
- ⌚ 37°C Temperature Controlled Assay Performance.
- ⌚ Low Reagent requirement (50ul per Test).
- ⌚ Measurement Starts when reagent is dispensed and Stops measurement when clot is detected. No need to press Start/Stop key manually.
- ⌚ Time, INR, Ratio displayed as a result.



Technical Specifications

Measuring System	: Photometric
Light Source	: LED
Display	: 4 Line x 16 Character, Blue LCD
Analysis	: PT, APTT
Results	: Time in seconds, INR, Ratio
Electrical Rating	: 12 VDC / 2 Amps

Operating Conditions

Temperature	: 5°C to 35°C
Relative Humidity	: Up to 85%
Size	: L:208 mm x W:108mm x H:43mm
Weight	: 0.500 Kg
Storage	: 10°C to 60°C, away from direct sunlight.
Usage	: Indoor use only

Due to continuous development, specifications are subject to change without prior notice.



STAR DELTA

641, Bhadrakali Compound,
Off Mumbai - Goa Highway,
A/P Zarap, Tal.: Kudal, Dist.: Sindhudurg,
Maharashtra - 416 520, INDIA
Tel : +91 2362 232030,
Email : stardelta216@gmail.com