



LABOMED®

ideas for vision

VISION 2000 Educational Microscope

Item

Specifications

Stand	Stable stand with convenient location of focus controls with special plastic hand rest
Viewing Body	Binocular tube , 45° inclined, 360° rotatable
Eyepiece	Wide field eyepiece 10x/18mm with foldable eye guard, Anti Fungus.
Nosepiece	Quadruple Nosepiece on ball bearing guideways with rubber grip.
Objectives	EP series DIN Semi Plan Achromatic objectives 4x, 10x, 40x (spring loaded), 100x (spring loaded, oil), anti fungus
Stage	Mechanical stage with ball bearing slides and 'soft feel' low drive co-axial control, with Graduated scale.
Condenser	Sub Stage abbe condenser NA 1.25 with aspheric lens. Iris diaphragm with snap in blue filter. Rack and Pinion movement on metal guides.
Focusing	Co-axial coarse and fine focusing mechanism system on ball bearing guide ways provided with Tension adjustment ring and auto focusing stop.
Illumination	Halogen 6V-20W with variable electronic regulator. Up to 2,000 hours of Halogen lamp life.
Electrical	Input 220V - 240V AC, CE Approved.
Accessories	Dust Cover, Day Light Mirror.

Configurations

Description	Configuration
Halogen	
Standard, 20W	Binocular



N. K. Jain Instruments Pvt. Ltd.
15, Industrial Estate, Ambala Cantt - 133 006
Haryana (INDIA)
info@labomedindia.com
www.labomedindia.com
Ph.: 91-171-2699015, 2699315
Fax: 91-171-2699316



Distributor